

Store Forward dermatology telemedicine is a convenient way to practice asynchronous consultation. Practitioners use their own computer system and must meet the following computer and network requirements to participate.

## Hardware Requirements

- Processor – Pentium 4 or better
- RAM – minimum 512MB of RAM
- Monitor – Minimum 1024 x 768 Resolution Monitor at 32 bit color depth
- Hard disk space – Minimum 100MB
- Free USB port – at least USB 2.0 (this is requirement for referring users only)

## Platform and Browser Requirements

- Windows 7, 8, 10 – Internet Explorer or Chrome
- Mac – Chrome or Safari
- iPhone and iPad – Safari
- Android – Chrome

Notes: (1) The *Medweb “Advanced Viewer”* is limited to Windows PC with Internet Explorer 11.

(2) The *OTN Telederm app* is available only to users with OTN credentials (the app is not integrated with ONE ID).

## Internet connectivity:

- Internet connectivity with at least 256 kbps upload bandwidth
- Allow secure access to Store Forward application server (*i.e.*, https://)

## Camera Requirements (referring users only):

- 8 megapixels or better
- Auto focus
- Macro mode (close up mode)
- Auto white balance

